



Anti- Non structural protein 1 (NS1) from Dengue virus 1+2+3+4 serotypes. (Llama serum)

Product code: ADVNS103

Description	Llama (<i>Llama glama</i>) anti-Dengue virus NS1 serum.
Specificity	NS1 from 4 Dengue virus serotypes (DV1+DV2+DV3+DV4).
Product Type	IgG polyclonal serum
Isotype, clonality	IgG Polyclonal
Conjugate	unconjugated
Purity	Whole serum
Form	Freeze dried.
Storage Buffer	N/A
Concentration	N/A
Size	100 uL
Storage instructions	Shipped with ice packs. Store at -20°C. Stable for several years.
Intended to use	RUO only, ELISA, WB.