



## Anti- Non structural protein 1 (NS1) from Dengue virus 1+2+3+4 serotypes. (Llama IgG)

*Product code: ADVNS104*

<b>Description</b>	Llama ( <i>Llama glama</i> ) anti-Dengue virus NS1 polyclonal IgG.
<b>Specificity</b>	NS1 from 4 Dengue virus serotypes (DV1+DV2+DV3+DV4)
<b>Product Type</b>	IgG polyclonal Antibody
<b>Isotype, clonality</b>	IgG, polyclonal
<b>Conjugate</b>	unconjugated
<b>Purity</b>	Affinity-purified antibody
<b>Form</b>	Freeze dried.
<b>Storage Buffer</b>	PBS 0,5M, Proclin
<b>Concentration</b>	1 mg/mL
<b>Size</b>	100 uL
<b>Storage instructions</b>	Shipped with ice packs. Store at -20°C. Stable for several years.
<b>Intended to use</b>	RUO only, ELISA, WB.