

## **Anti-native Flagellin from *C. chauvoei* biotinylated V<sub>H</sub>H domain recombinant antibody**

*Product code: HCMRA10*

<b>Description</b>	Biotinylated <i>Clostridium chauvoei</i> FliC recombinant binding protein (V <sub>H</sub> H domain)
<b>Specificity</b>	<i>Clostridium chauvoei</i> FliC protein
<b>Product Type</b>	Biotinylated Recombinant single-domain antibodies (V <sub>H</sub> H) produced in <i>E. coli</i> .
<b>Isotype</b>	V <sub>H</sub> H domain of alpaca heavy-chain antibody, monoclonal, monovalent
<b>Conjugate</b>	Biotin (monobiotinylated)
<b>Purity</b>	Affinity-purified antibody fragment
<b>Form</b>	Freeze dried.
<b>Storage Buffer</b>	PBS 0,5M, BSA 1%. Proclin.
<b>Concentration</b>	1mg/mL
<b>Size</b>	0.1mg
<b>Storage instructions</b>	Shipped with ice packs. Store at -20° C. Stable for several years.
<b>Intended to use</b>	RUO only, ELISA, EIA.