



Anti-native flagellin from *C. chauvoei* hyperimmune serum

Product code: AFC00

Description	Rabbit anti- <i>C. chauvoei</i> flagella hyperimmune serum
Specificity	Flagellin protein from <i>C. chauvoei</i> IRP 631 (USDA-APHIS).
Product Type	Hyperimmune serum
Isotype	Undetermined mainly IgG, Polyclonal
Conjugate	unconjugated
Purity	Bulk serum
Form	Bulk, freezed.
Storage Buffer	NA
Concentration	ND
Size	100 µL
Storage instructions	Shipped with ice packs. Store at -20°. Stable for several years.
Intended to use	RUO only, ELISA, EIA, WB.