

## Anti-R. equi VapA V<sub>H</sub>H domain recombinant antibody

Product code: HCMRA14

<b>Description</b>	<i>Rhodococcus equi</i> Virulence Associated Protein A (VapA) recombinant binding protein (V <sub>H</sub> H domain)
<b>Specificity</b>	R. equi VapA protein
<b>Product Type</b>	Recombinant single-domain antibodies (V <sub>H</sub> H) produced in E. coli.
<b>Isotype</b>	V <sub>H</sub> H domain of alpaca heavy-chain antibody, monoclonal, monovalent
<b>Conjugate</b>	unconjugated
<b>Purity</b>	Affinity-purified antibody fragment
<b>Form</b>	Freeze dried
<b>Storage Buffer</b>	PBS 0,5M, BSA 1%. Proclin
<b>Concentration</b>	1mg/mL
<b>Size</b>	20,1mg
<b>Storage instructions</b>	Shipped with ice packs. Store at -20°. Stable for several years.
<b>Intended to use</b>	RUO only, ELISA, EIA..