

## Anti-canine LH V<sub>H</sub>H domain recombinant antibody

*Product code: HCMRA2*

<b>Description</b>	Canine Luteinizing Hormone recombinant binding protein (V <sub>H</sub> H domain)
<b>Specificity</b>	Canine Luteinizing Hormone
<b>Product Type</b>	Recombinant single-domain antibodies (V <sub>H</sub> H) produced in <i>E. coli</i> .
<b>Isotype</b>	V <sub>H</sub> H domain of alpaca heavy-chain antibody, monoclonal, monovalent
<b>Conjugate</b>	Unconjugated
<b>Purity</b>	Affinity-purified antibody fragment
<b>Form</b>	Freeze dried.
<b>Storage Buffer</b>	PBS 0,5mM
<b>Concentration</b>	1mg/ml
<b>Size</b>	0.1mg
<b>Storage instructions</b>	Shipped with ice packs. Store at -20°C. Stable for several years.
<b>Intended to use</b>	RUO only, ELISA, EIA, WB.