



Depleted Anti-native flagellin from *C. chauvoei* IgG (rabbit)

Product code: AFC02

Description	Rabbit anti- <i>C. chauvoei</i> FliC polyclonal IgG depleted from IgG against <i>C. septicum</i> FliC.
Specificity	Anti- FliC flagellin protein from <i>C. chauvoei</i> IRP 631(USDA APHIS strain)
Product Type	IgG polyclonal Antibody Purified using Protein A
Isotype	IgG polyclonal.
Conjugate	unconjugated
Purity	Affinity-purified antibody fragment
Form	Freeze dried
Storage Buffer	PBS 0,5M, Proclin.
Concentration	1mg/mL
Size	0,1 mL
Storage instructions	Shipped with ice packs. Store at -20°. Stable for several years.
Intended to use	RUO only, ELISA, EIA, WB.