



Anti-recombinant NS1 from Zika virus (option 2) Biotinylated V_HH domain recombinant antibody

Product code: HCMRA12

Description	Biotinylated Zika Virus, Nonstructural Protein 1 (ZIKV, NS1) recombinant binding protein (V _H H domain) (option 2)
Specificity	Zika virus NS1 protein
Product Type	Biotinylated Recombinant single-domain antibodies (V _H H) produced in E. coli. His-Tag, HA-Tag.
Isotype	V _H H domain of alpaca heavy-chain antibody, monoclonal, monovalent
Conjugate	Biotin (monobiotinylated)
Purity	Affinity-purified antibody fragment
Form	Freeze dried.
Storage Buffer	PBS 0,5M, BSA 1%. Proclin.
Concentration	1mg/mL
Size	0.1mg
Storage instructions	Shipped with ice packs. Store at -20° C. Stable for several years.
Intended to use	RUO only, ELISA, EIA.