



Anti- Non structural protein 1 (NS1) from Dengue virus 1+2+3+4 serotypes. (Rabbit IgG)

Product code: ADVNS102

Description	Rabbit (<i>Oryctolagus cuniculus</i>) anti-Dengue virus NS1 polyclonal IgG.
Specificity	NS1 from 4 Dengue virus serotypes (DV1+DV2+DV3+DV4)
Product Type	IgG polyclonal Antibody
Isotype, clonality	IgG, polyclonal
Conjugate	unconjugated
Purity	Affinity-purified antibody
Form	Freeze dried.
Storage Buffer	PBS 0,5M, Proclin
Concentration	1 mg/mL
Size	100 uL
Storage instructions	Shipped with ice packs. Store at -20°C. Stable for several years.
Intended to use	RUO only, ELISA, WB.