



Anti-BSA V_HH domain biotinylated recombinant antibody

Product code: HCMRA9

Description	Biotinylated Bovine Serum Albumin recombinant binding protein (V _H H domain)
Specificity	Bovine Serum Albumin (BSA)
Product Type	Biotinylated Recombinant single-domain antibodies (V _H H) produced in E. coli. His-Tag, HA-Tag.
Isotype	V _H H domain of Llama heavy-chain antibody, monoclonal, monovalent
Conjugate	Biotin (monobiotinylated)
Purity	Affinity-purified antibody fragment
Form	Liquid solution.
Storage Buffer	PBS 0.5M, Proclin 0.1%
Concentration	1mg/mL
Size	0,1mg
Storage instructions	Shipped with ice packs. Store at -20°. Stable for several years.
Intended to use	RUO only, ELISA, EIA.