



Anti-canine GnRH biotinylated V_HH domain recombinant antibody

Product code: HCMRA7

Description	Biotinylated Canine GnRH recombinant binding protein (V _H H domain)
Specificity	Canine Gonadotropin Realizing hormone (GnRH).
Product Type	Biotinylated Recombinant single-domain antibodies (V _H H) produced in E. coli. His-Tag, HA-Tag.
Isotype	V _H H domain of Llama heavy-chain antibody. monoclonal, monovalent
Conjugate	Biotin (monobiotinylated)
Purity	Affinity-purified antibody fragment
Form	Freeze dried.
Storage Buffer	PBS 0,5M, BSA 1%. Proclin.
Concentration	1mg/ml
Size	0.1mg
Storage instructions	Shipped with ice packs. Store at -20°C. Stable for several years.
Intended to use	RUO only, ELISA, EIA.