



## Anti-Progesterone polyclonal antibody

*Product code: POLY01*

<b>Description</b>	Rabbit ( <i>Oryctolagus cuniculus</i> ) anti-Progesterone polyclonal IgG.
<b>Immunogen</b>	Progesterone-3-CMO-BSA Antigen
<b>Product Type</b>	IgG polyclonal Antibody Purified using Protein A
<b>Isotype</b>	IgG Polyclonal
<b>Conjugate</b>	unconjugated
<b>Purity</b>	Affinity-purified antibody fragment
<b>Form</b>	Freeze dried.
<b>Storage Buffer</b>	PBS 0,5M, Proclin
<b>Concentration</b>	1mg/ml
<b>Size</b>	100uL
<b>Storage instructions</b>	Shipped with ice packs. Store at -20°. Stable for several years.
<b>Intended to use</b>	RUO only, ELISA.