



Anti- Non structural protein 1 (NS1) from Dengue virus 1+2+3+4 serotypes. (Rabbit serum)

Product code: ADVNS101

| | |
|-----------------------------|--|
| Description | Rabbit (<i>Oryctolagus cuniculus</i>) anti-Dengue virus NS1 serum. |
| Specificity | NS1 from 4 Dengue virus serotypes (DV1+DV2+DV3+DV4). |
| Product Type | IgG polyclonal serum |
| Isotype, clonality | IgG Polyclonal |
| Conjugate | unconjugated |
| Purity | Whole serum |
| Form | Freeze dried. |
| Storage Buffer | N/A |
| Concentration | N/A |
| Size | 100 uL |
| Storage instructions | Shipped with ice packs. Store at -20°C. Stable for several years. |
| Intended to use | RUO only, ELISA, WB. |